

HOME CONTAINER
Daily Care Instructions



Do not open the chamber or handle the proxy directly.

All maintenance is performed through sealed cartridges and automated monitoring.

TARGET CONDITIONS


37.0°C

5.0% CO₂

94% RH

pH 7.4

CO₂ SUPPLY
Provides a stable 5.0% CO₂ environment.



1 PLACEMENT

- Place on a stable, level surface.
- Keep at least 20 cm clearance around the device.
- Reconcentrate room temperature: 18–26°C.
- Do not move the container while a institution or circulation cycle is active.



2 POWER & CONNECTIVITY

- Keep the container connected to power.
- Use only the supplied power adapter.
- Connect to Wi-Fi for monitoring, software updates, and clinic alerts.
- The internal backup battery supports the proxy for approximately 60 minutes.



3 NUTRITION

- Check the medium level once daily.
- Replace the cartridge when the screen displays REPLACE medium.
- Do not refill, dilute or reuse a cartridge.
- Scan the new cartridge before installation.



Used cartridges must be sealed in the supplied return pouch and returned to the clinic.



4 LIGHTING

- Keep away from windows, and direct sunlight.
- Maintain a 12 h light / 12 h dark cycle.
- Avoid direct spotlight and illumination.



5 DAILY TRAINING

- Complete 3 guided sessions per day
- Morning session Afternoon session Evening session
- Each session lasts 8–12 minutes.
 - leave at least three hours between sessions.




6 DAILY CHECK

- Confirm the following twice a day:
- power is active
 - Network connection stable
 - Temperature, humidity, CO₂ and pH are within range.
 - Medium level is sufficient
 - No tubing is bent, disconnected or leaking
 - No warning message.

7 TERMINATION

Termination is possible at any stage following clinic authorization



Press and hold **TERMINATION** on the device interface for 10 seconds.

Wait until the screen confirms **BIOLOGICAL ACTIVITY ENDED**.

Do not open the chamber or remove the proxy.

Disconnect the device only after the system displays **SAFE TO TRANSPORT**.

Return the sealed container to clinic.



DO NOT DISPOSE OF THE PROXY, NUTRIENT MEDIUM, CARTRIDGES, OR TUBING IN HOUSEHOLD WASTE



CLINICAL SUPPORT
24/7



+972-3-1234567



support@proxybond.clinic



proxybond.clinic/support